



SIMPLIFY
THE POLISHING
PROCESS



MAGNETIC BACKED Polishing Cloths





Save Time

Save valuable time by eliminating the process of removing a worn polishing cloth and adhering a new cloth to a carrier plate. Magnetic backed polishing cloths are ready to use straight out of the box.

Repeatable Results

Ensure repeatability in your polishing process. Buehler cloths use a stiff and durable magnetic backing to provide repeatable polishing results over the life of the cloth.




Advanced Solutions

Buehler's offering of polishing cloths provides industry leading solutions for many different applications. Trust the experts at Buehler to develop a process to meet the specific needs of your lab.

Available on Your Favorite Polishing Cloths

A magnetic receptive backing is available on most of Buehler's polishing cloth offering. These products combine the original Buehler polishing cloth material with a stiff magnetic liner to provide a convenient and easy to handle solution.



Ordering Information

Cloth	8in [203mm]	10in [254mm]	12in [305mm]
UltraPad	42-3008	42-3010	42-3012
Nylon	42-3108	42-3110	42-3112
TexMet C	42-3208	42-3210	42-3212
TriDent	42-3308	42-3310	42-3312
VerduTex	42-3408	42-3410	42-3412
VelTex	42-3508	42-3510	42-3512
PoliCloth	42-3608	42-3610	42-3612
MicroCloth	42-3708	42-3710	42-3712
MicroFloc	42-3808	42-3810	42-3812
MasterTex	42-3908	42-3910	42-3912
ChemoMet	42-4008	42-4010	42-4012

Each item number above contains 5 pcs



BUEHLER Worldwide Headquarters
North America-South America Offices
 41 Waukegan Road
 Lake Bluff, Illinois 60044-1699 USA
 P: 800 BUEHLER (800-283-4537)
 P: 847-295-6500
www.buehler.com | info@buehler.com

European Headquarters
BUEHLER Germany
info.eu@buehler.com

BUEHLER France
info.eu@buehler.com

BUEHLER United Kingdom
info.eu@buehler.com

BUEHLER China
info.cn@buehler.com

BUEHLER Japan
info.japan@buehler.com

BUEHLER Asia-Pacific
info.asia@buehler.com