



# EpoKwick™ FC

## Fast Curing Premium Epoxy System

Spend **less** time preparing,  
and **more** time analyzing.

### Efficient

No need to wait for the highest quality mounts. Cure your mounts without the use of an oven in only 2 hours.



### Maximum Pore Penetration

EpoKwick FC has the lowest mix viscosity available in a room temperature cured epoxy. Lower viscosity improves sample quality by maximizing penetration of gaps, cracks and pores.

### Excellent Edge Retention

EpoKwick FC has strong adherence and zero shrinkage with most samples. This means no gaps between the mount material and your sample, protecting the sample edge and minimizing contamination.

# Obtain The Highest Sample Preparation Quality.


EpoKwick FC is a unique clear epoxy system formulated specifically for applications needing high quality mounts in the fastest time. EpoKwick FC has the lowest mix viscosity and the lowest shrinkage on the market combined with one of the fastest cure times available without the use of an oven.

	<b>EpoKwick™ FC</b>	<b>Nearest Competitor*</b>
<b>Viscosity</b>	250 cP	430 cP
<b>Shrinkage</b>	<1 um	5 um
<b>Cure Time</b>	2 hours	12 hours
<b>Peak Exotherm</b>	250° F (121° C)	115° F (46° C)
<b>Shore D Hardness</b>	82	78

\*Comparison of product under the same standard test conditions. Results may differ from competitor's published specifications.



Inside diameter of a steel washer mounted with EpoKwick FC showing no shrinkage gap.



Inside diameter of a steel washer mounted with Nearest Competitor epoxy produces shrinkage gaps.

## Ordering Information

Part Number	Description
20-3453-128	128 oz (3.8L) Resin
20-3453-032	32 oz (0.95L) Hardener



While EpoKwick FC is excellent for most castable mounting applications, it is not suitable for temperature sensitive samples such as plastics. If your sample is temperature sensitive use Buehler EpoThin™ 2 or EpoxiCure™ 2.



**BUEHLER Worldwide Headquarters  
North America-South America Offices**  
41 Waukegan Road  
Lake Bluff, Illinois 60044-1699 USA  
P: 800 BUEHLER (800-283-4537)  
P: (847) 295-6500  
www.buehler.com | info@buehler.com

**European Headquarters  
BUEHLER Germany**  
info.eu@buehler.com

**BUEHLER France**  
info.eu@buehler.com

**BUEHLER United Kingdom**  
info.eu@buehler.com

**BUEHLER China**  
info.cn@buehler.com

**BUEHLER Japan**  
info.japan@buehler.com

**BUEHLER Asia-Pacific**  
info.asia@buehler.com