

MICROSCÓPIO PORTÁTIL MODELO BXJ 2000.



BXJ2000

This model BXJ2000 portable metallurgical microscope, it is comfortable to use and small volume, it can be with the metallurgical sample tool.. Total magnification can be from 100X to 500X.

The C-Mount like the magnification 0.5X/0.7X/0.35X are optional.

Specification	BXJ2000		
Eyepiece	WF10X/18mm		
	WF12.5X/16mm for optional		
Plan Achromatic Objective	40X	20X	10X
Magnification	400X	200X	100X
NA	0.25	0.12	0.10
Working Distance(mm)	7	18	18
Stage with Magnet	Stage size: 100×80mm, moving range: 10×10mm		
Focusing System	Separate coarse and fine. Coarse range 30mm and fine range: 1mm		
Illumination	Halogen lamp 12V 20W		